



# HCHS/SOL Anthropometry Quality Control

ID NUMBER:							
------------	--	--	--	--	--	--	--

FORM CODE: AQC  
VERSION: 2, 02/14/2022

Contact  
Occasion

0	3
---	---

Occurrence

0	1
---	---

## Administrative Information

0a. Completion Date: / /

0b. Staff ID:

0c. Anthropometry QC selection  (0=NOT selected, 1= Selected)

**[If participant is NOT selected for AQC (Q0c=0), SAVE and close the form]**

0d. History of implanted electronic heart pacemaker or defibrillator (0=No, 1=Yes)  **[Prefill from PSE3]**

**[If 0d=1 you MUST set Tanita scale to weigh-only mode to weigh participant]**

**Instructions:** Enter the answer given by the participant for each response. Set CDART Field Status to 'Refused', 'No Response', 'Missing' etc. for those questions that do not list these values as possible answer choices. In order to measure bioimpedance, the participant must be barefoot. Set the Tanita analyzer to report metric units (cm/kg).

## A. DETERMINATION OF ABILITY TO STAND

1. Assessment of ability to stand (choose one):
- Can stand erectly on both feet. .... 1
- Can stand on both feet, but posture not erect. .... 2
- Cannot stand on both feet. .... 3 → **GO TO ITEM 10**

## B. HEIGHT, WEIGHT, and BIO-IMPEDANCE

2. Standing height (round to nearest cm):    cm
3. a) Self-reported weight (to the nearest lb or kg):
- b) Units (check one):  lb  kg
4. Weight:    .  kg
5. Fat (%):   .  %
6. Impedance:    Ohms
7. Fat mass:    .  kg
8. Lean body mass (FFM):    .  kg
9. Total body water (TBW):    .  kg

## C. BODY SIZE

10. Girth (round to nearest cm)
- a) Waist:    cm

ID NUMBER:							
---------------	--	--	--	--	--	--	--

FORM CODE: ANT  
VERSION: A  
7/06/07

Contact Occasion		
---------------------	--	--

SEQ #		
-------	--	--

b) Hip: .....

--	--	--

 cm